

Revolutionary MicroSpring™ textile offers breathability, moisture transfer and anti-microbial protection not possible with gels, foams or fleece.

Microfilaments flex to reduce pressure on bony prominences such as the heel, elbow, vertebra and trochanter.

Layers of MicroSpring Textile puts up to 20 thousand microfilaments per square inch to act as mini pressure relieving springs.



- Superior skin pressure and shear relief
- Breathable - fabric and filament layers allow air flow
- Moisture wicking - allows moisture transfer
- Can be washed and tumble dried up to 95°C (200°F) for infection control and allows for multi-patient use
- Durable - will not flatten out or become flaky
- Stretchy knit (not woven) fabric has "give" to reduce surface friction
- Anti-microbial - will not support growth of bacteria
- Compatible with topical creams and ointments
- Latex free - 92% polyester, 8% spandex
- CA117 Flame Resistant

Durometer reading from skin (30)

Durometer reading from DermaSaver (30) matches skin



Unlike foams, gels, fleece and polyfil, DermaSaver™ reduces all of the conditions that cause skin breakdown.

Breathability And Moisture Transfer
Ample ventilation space between the microfilaments and textile layers allow for air flow and moisture transfer.

Microspring Textile
Microfilaments flex to reduce pressure on bony prominences such as the heel, elbow, vertebra and trochanter. Layers of Microspring® Textile put thousands of microfilaments per square inch to act as mini pressure reducing springs.

Protected by the ÆGIS Microbe Shield
All layers of DermaSaver products are protected by the aegis microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product. It is permanently bonded to the fibres, does not leach, or volatilize.



Totally Launderable
Machine wash & dry even in the high temperature of institution laundries. Allows for multi patient use. Lower wash/dry temperatures will lengthen product life. If using bleach, use lower concentrations.

“ We’ve had problems with other protectors sticking to the wounds, but not with DermaSaver. I really like the air permeability and anti-microbial qualities of the products, too. ”


- Diane Husar, Dir. Quality Assurance, United Methodist Homes.



HealthSaver Pty Ltd,
14/140 Wecker Road, Mansfield, QLD 4122, Australia.
Phone: 1300 767 888 Fax: 1300 767 999
Email: info@healthsaver.com.au
www.healthsaver.com.au



Unlike foams, gels, fleece and polyfil DermaSaver™ reduces ALL of the conditions that cause skin breakdown. Protect all these areas that are vulnerable to pressure sores.



Breathable • Anti-microbial • Durable • Washable

Available exclusively from HealthSaver
1300 767 888
www.healthsaver.com.au



DermaSaver SkinTubes™ for Arms & Legs

- Superior Protection From Skin Tears
- Non-Constricting
- Anti-Microbial
- Breathable
- Totally Launderable



These soft flexible tubes fit loosely and remain elongated on arms & legs.

Two layers of soft, lightweight MicroSpring Textile™ form a loose-fitting tube that stays in place to protect fragile skin from pressure, friction and minor traumas that result in tears, bruising and minor abrasions. Unlike traditional fabric stockings, which are tight and binding, the layers of microfilaments allow air-flow and moisture transfer.

Simply slip the tube over the arm or leg, there are no hook and loop fasteners or tight fitting spandex. The narrower end of the tube gently rests on the ankle or wrist to keep it from falling off. The cross microfilaments in the MicroSpring Textile allows these soft tubes to be self-supporting, which means they won't slip down or bunch up.

SkinTubes are protected by the Aegis Microbe shield which uses an EPA registered anti-microbial technology to inhibit the growth of micro-organisms on the product.

You can use creams and ointments with SkinTubes. All SkinTubes can be completely washed and dried even in high temperature institution laundries.

SkinTubes are also available with double protection for knees or elbows.

DermaSaver SkinTubes™ for Legs



DermaSaver SkinTubes™ for Arms



The DermaSaver Arm & Leg tubes are great products for skin protection. Residents find them soft and comfortable, and not tight like other products we've tried.

- Loretta Gerving, DON. Marian Manor Healthcare Centre.

DermaSaver Heel & Ankle Protection



DermaSaver Speciality Items

