

ETON

'Box in box' dual motor lift chair with headrest



* Shown in Stonewash - Amber

Description and Features



Aside from its stylish appearance, the Eton also boasts the new integral 'box in box' mechanism. This fully enclosed feature means that nothing can inadvertently obstruct the chair in its fully lifted position. The Eton also includes an independent headrest function for enhanced comfort.

- New 'box in box' technology meaning a greater capacity
- Two motor lift recliner
- Independent headrest function
- Handwand with USB port

Specifications

Weight Capacity	226.7 kgs (500 lbs)
Seat Width	53.3 cm (21")
Seat Depth	52.0 cm (20.5")
Chair Width	87.6 cm (34.5")
Chair Depth	95.3 cm (37.5")
Chair Height	106.7 cm (42")
Operation	Handwand

Fabric Selection

-  Stonewash - Amber
-  Stonewash - Dune