

FUSION CHAIR

The Fusion is a high back chair with multi-layered seat and back, shaped for comfort and support.



FEATURES

- The Fusion was developed specifically with high care clients in mind.
- The anti-bacterial cover is made from an easy care water resistant 4-way stretch fabric, suitable for sufferers of incontinence.
- Supportive Dunlop Visco elastic low resilience foam contours perfectly to the body shape helping reduce pressure points.
- The ultimate in comfort, support and adjustability with a range of custom sizes available.
- Available in Renaissance Vinyl and Fusion.

SPECIFICATIONS

CODE	SEAT WIDTH	MIN. SEAT DEPTH	MAX. SEAT DEPTH	MIN. SEAT HEIGHT	MAX. SEAT HEIGHT	BACKREST HEIGHT	MAX. USER WEIGHT
18030	520mm	500mm	650mm	450mm	550mm	1130mm	250KG
18030B65	650mm	500mm	650mm	450mm	550mm	1130mm	250KG
18030BC	Custom	500mm	650mm	500mm	500mm	1180mm	250KG

R&R Healthcare Equipment

FRECALL 1800 684 405 Facsimile 03 5021 5653

Email sales@randrhealthcare.com.au

www.randrhealthcare.com.au



Innovative Healthcare Solutions