

COMMODE SEAT KA11HOP-B

NEW fabric and colour range

FEATURES

- Anti-fungal and anti-bacterial material
- Easily maintained and wiped clean
- Superior tear strength with great flexibility
- Waterproof
- Abrasion tested to 100,000+ cycles
- Clip-on design including pan carrier runners is easily removed for cleaning and service.
- Interchangeable with other models to suit individual occupants.
- Can be fitted to all 445mm wide commodes.
- White vinyl available as special order



Fabric



Anti fungal



Fire resistant



Anti bacterial



Wear resistant tested

SPECIFICATIONS

Description	Overall Width	Overall Depth	Overall Height	Aperture Width	Aperture Depth	Product Weight	ARTG
Robin Blue Vinyl Seat KA11HOP-B	445mm	430mm	55mm	200mm	260mm	2kg	209826

Improving your
life everyday™

kcare.com.au

KCare
Healthcare Solutions