

# Shoprider COUGAR

ASK FOR THE BRAND SHOPRIDER



**COUGAR 14**



**COUGAR 10**

**SHOPRIDER®**  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## SHOPRIDER<sup>®</sup> mobility COUGAR

www.shoprider.com.au

The SHOPRIDER™ COUGAR is ergonomically engineered to provide the ultimate combination of style and comfort. **Cougar** with its powerful smooth motors, mid wheel maneuverability, and comfortable seating provides a safe and comfortable ride. The smooth operation combined with mid wheel drive means easy operation under various conditions. It's ideal for both inside and outside use on even ground and can travel long distances between charges.

SPECIFICATIONS	G424M	G424L	FEATURES:
Length	920mm	920mm	Large Captain's Seat
Width	620mm	690mm	Retractable charger cord
Height	1120mm	1180mm	Adjustable arm rests
Battery Capacity (Minimum Amp)	2 x 40Amp	2 x 40Amp	Option of swing away leg rest
#Maximum Range	up to 30km	up to 40km	Option elevating leg rest
Maximum Load Single User*	136Kgs	136Kgs	Sleek body design
Charger on board	4Amp	8Amp	Color change panels
Motor Size	0.8 Hp	1.34 Hp	Smooth operation
Maximum Speed	8Kph	9.6Kph	Free wheel device
*Safe Climbing Angle	8 Degrees	8 Degrees	Variable seat height
*Max Climbing Angle	12 Degrees	12 Degrees	Controller interchangeable both arms
Turning Radius	570mm	600mm	Rear drive bracket available
Wheel Size:			
- Front Caster	2 x 125mm	2 x 125mm	
- Middle (drive)	2 x 250mm	2 x 350mm	
- Rear caster	2 x 100mm	2 x 100mm	


Driver weight may exceed weight of unit. Reduce speed when turning

### IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



**SHOPRIDER<sup>®</sup>**  
www.shoprider.com.au  
mobility

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au

Your Authorised dealer

**LAKESIDE**  
mobility

2/1 Metier Linkway | Birtinya | QLD | 4575  
(07) 5325-1231 (07) 5391-4012