

Shoprider

PUMA

ASK FOR THE BRAND SHOPRIDER



PUMA 10



PUMA 14HD

SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPPRIDER®
www.shoprider.com.au *mobility*

PUMA

The SHOPPRIDER™ PUMA has been ergonomically engineered to provide a robust combination of style and comfort. The PUMA is fitted with powerful smooth motors, mid wheel manoeuvrability, and comfortable seating. The smooth operation of the front and rear elastometers combined with mid wheel drive means easy operation on a variety of terrains and hard surfaces. It's ideal for both inside and outside use on even ground and can travel long distances between recharging.

SPECIFICATIONS	PUMA10	PUMA14HD	FEATURES
Length	1060mm	1060mm	Large Captain's Seat
Width	620mm	640mm	20" Seat PUMA 14HD
Height	1160mm	1180mm	Sleek body design
Battery Capacity	2 x 36Amp	2 x 75Amp	Smooth operation
#Maximum Range	up to 20km	up to 25km	Variable Seat height
Maximum Load Single User*	136Kgs	205Kgs	Controller interchangeable both arms
Charger on board	4Amp	8 Amp	Free wheel device
Motor Size	0.8 Hp	1.34 Hp	Adjustable arm rests
Maximum Speed	8Kph	6Kph	Option elevating leg rests
*Safe Climbing Angle	8 Degrees	8 Degrees	Option of swing away leg rest
*Max Climbing Angle	12 Degrees	9 Degrees	Stamp supports
Turning Radius	570mm	600mm	Comprehensive Warranty
Wheel Size:			
- Front Caster	2 x 175mm	2 x 175mm	
- Middle (drive)	2 x 260mm	2 x 350mm	
- Rear caster	2 x 175mm	2 x 175mm	


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPPRIDER®
www.shoprider.com.au *mobility*

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au

Your Authorised dealer

LAKESIDE
mobility

2/1 Metier Linkway | Birtinya | QLD | 4575
(07) 5325-1231 (07) 5391-4012