

NEW!

GIVE YOURSELF A NEW LIFE, INDOORS AND OUTDOORS.

The Let'sGo Indoor Walker has become a great success around the world. An indoor walker meeting all the demands regarding quality and environment, support and safety. ISO certified. Easy to use with numerous clever features. A helpful companion at home, at the hospital or in retirement homes.



MORE ACCIDENTS OCCUR AT HOME

Most accidents among the elderly happen at home and the number of accidents increases rapidly with older age. Falling accidents at home can result in serious injury.

Why an Indoor Walker?

At old age many need walkers to be able to move safely outdoors, but few want to use them indoors as they are big, clumsy, impractical and bring dirt into the home. So a walking stick is often used which does not replace the safety that a walker provides.

Now there is a practical solution which has been successful around the world; the Let'sGo Indoor Walker, created in Sweden and fulfilling requirements around the world. It has passed all the tests for quality and environment.

The Benefits!

It is small and easy to use in narrow spaces, turns easily and gives support and safety in all situations.

Width only 54cms and weighs only 6.5kgs.

The handle is adjustable from 79 to 90 cms to suit different heights.

Easy to manage brakes that may be operated with one hand only and the one brake operates both wheels. And of course a parking brake.

With a removable plastic tray, use it as a help in the kitchen as well as for transporting food & drinks or as a small table to eat meals from.

Use it as a support if you need to get up from bed in the middle of the night. Take it to the bathroom & toilet for added safety.

It is beautiful and elegant in all detail, but the most important thing is that it helps many to a new and more comfortable life.

Worldwide Success.

The Let'sGo has been a great success internationally. Today it is sold in 17 countries from Germany to Australia. Selling over 20,000 units in its 1st Year.

Beech 0020

Red & Black 0030



Let'sGo Indoor Walker

Measurements & Weight

Handle	79-90cm
Tray Height	68cm
Length	63cm
Width	54cm
Weight	6.5kgs

Practical Fabric Basket Let'sGo

Accessory basket that is fixed under the tray with a Velcro band. Max weight in basket 3kgs



Folding

Easy to fold and transport.



DVA Approved*

TRUST CARE



New in 2012!

Let'sGo Out Gives Quality of Life.

Quality of life is to be able to walk outdoors under safe and comfortable conditions, simply and easily. Therefore we have developed the Let'sGo Out further to the Let'sGo Indoor Walker. Many of the features are the same but a few extra features have been added.

Let'sGo Out is Filled with More Features.

It is light and elegant, both folded and unfolded. It has an X frame cross folding system which makes it easy to fold, takes up little space and stands alone when folded.

Big front and rear wheels make it easy to get over curbs and other obstacles. Even more important is the spring system in the wheels that gives a smooth movement compensating for stiff joints.


The Let'sGo Out has soft, comfortable rubber handles that are adjustable in height from 78 to 96 cms to fit different heights. It has an easy to manage brake which of course can be locked and a comfortable seat in imitation leather for strength.

The Let'sGo Out has an elegant design in every aspect.

It will come in two colors: 0120 Black/Silver/Black Seat.

0150 Bronze/Silver/Brown Seat.

-Bag & Back-strap included 0140 Tray Accessory



Measurements & Weight

Handle	78-96cm
Handle Width	41cm
Seat Height	59cm
Length	64cm
Width	57cm
Weight	6.3kgs

Folding


Easy to fold and transport,
Stands by itself when folded.

Let'sGo Out Rollator

Tested and Approved

The Walkers are tested
according to:

ISO 11199-2:2006, approved
as a walker for people
weighing up to 130kgs.



Also CE-marked. 

Your Retailer is:

LAKE SIDE
mobility

2/1 Metier Linkway | Birtinya | QLD | 4575
☎ (07) 5325-1231 📠 (07) 5391-4012