

# Peak CK835 STD Electric Homecare Bed

Product Code: 9346376012003

## Features:

- Powder Coated Steel Frame
- 4 Locking Castor Wheels
- Detachable Timber Head/Foot Board and Full Length Side Panel
- Collapsible Side Rails
- 3 Motors
- Tested for AS/NZS 3200.1.0 & AS/NZS 3200.2.38, Section 4, Hospital Beds



### Product Specifications

Length	205cm
Width	99cm
Height Range	41cm-75cm
Capacity	220kg
Product Weight	84kg

### Product Specifications

Backrest Adjustment	0 - 80°
Knee- Break Adjustment	0 - 4°
Wheel Size	12.7cm/5"

## ARTG Listing:

- 206378