



oscar
FURNITURE


Easy Glide Diner

Dining made easy for everyone



Features:

- Allows the user or carer to effortlessly roll into the table
- 50mm Rear self locking wheels (release when seated locks when you stand for stability)
- Light and Rigid Aluminium Frame 38mm front fixed wheels
- Easy clean wood grain finish
- 120kg weight capacity
- High resilience foam
- Oscar Dot fabrics available for fast delivery - other fabrics available on request
- Available with matching arm pads



Oscar Dot Range

Specifications:

Seat Height	Seat width	Seat Depth	Chair Height	Back Height	Chair Width	Arm Height
480mm	420mm	450mm	830mm	370mm	570mm	690mm