

THINK OF COMFORT

TOPFORM™

Est. 1965



**CHOOSE
YOUR
MATERIAL**

- FABRIC
- VINYL
- LEATHER

**CHOOSE
YOUR
FUNCTION**

- Manual Recliner
- Electric Recliner
- 1 Motor Lift Recliner
- 2 Motor Lift Recliner



Petite



Mini



Medium



Maxi



Tall

SIZE	PETITE	MINI	MEDIUM	MAXI	TALL
Seat Height	47cm	47cm	50cm	53cm	54cm
Seat Depth	45cm	49cm	50cm	52cm	53cm
Seat Width	44cm	48cm	48cm	53cm	54cm
Seat To Top Of Back	57cm	60cm	64cm	68cm	74cm
Chair Width	78cm	78cm	78cm	83cm	83cm
Chair Height	99cm	102cm	107cm	117cm	123cm



One Motor



Two Motor

1 MOTOR LIFT RECLINER

- Wall saver mechanism
- Elevates feet fully before reclining backrest
- Easy 2 button handset
- Waterfall backrest with zips for adjustability
- 155kg capacity

2 MOTOR LIFT RECLINER

- Independent footrest and backrest functions
- Lies fully flat
- Waterfall backrest with zips for adjustability
- 155kg capacity.

OPTIONAL
COVERS

OPTIONAL
JESCA
PILLOW

OPTIONAL
INCONTINENCE
COVER

OPTIONAL
FOOTREST
CUSHION

OPTIONAL
LOCKING
CASTERS

OPTIONAL
ROUND
DRINK
SERVER

Some customisation available. Please check with sales staff

Also Available as pressure
reducing chair Luxor



With or without removable cushion
insert in seat to fit pressure reducing
air cushion such as ROHO.

Your local distributor:



2/1 Metier Linkway | Birtinya | QLD | 4575
 ☎ (07) 5325-1231 📠 (07) 5391-4012