

Ellipse 6" Petite Rollator

Product Code: 8217-FX













Features:

- ▶ AU STD. Tested 491649-1
- ▶ PU Seat With No Sharp Edges
- ▶ Easy to Grip Brakes
- ▶ Foldable
- ▶ Long Wearing Finish



ARTG Listing:

- ▶ 186869

Length (cm)		60.5	Folded Width (cm)		20	Front Castor Diameter (cm)		15.2/6"
Width (cm)		60.5	Seat Depth (cm)		31	Rear Castor Diameter (cm)		15.2/6"
Max Height (cm)		93	Seat Width (cm)		35.8	Max Load (kg)		150
Minimum Height (cm)		72	Seat Height from the Ground (cm)		47	Net Weight (kg)		7