

Peak Taima Aluminium Rollator - Small

Product Code: C4505-A-FX













Features:

- ▶ One Hand Folding Mechanism
- ▶ Secure Brake Design
- ▶ Mesh Bag Included
- ▶ Cane Holder Included
- ▶ Foldable



ARTG Listing:

- ▶ 227601

Length (cm)	 67	Folded Width (cm)	 22	Front Castor Diameter (cm)	 20.3/8"
Width (cm)	 44	Seat Depth (cm)	 22	Rear Castor Diameter (cm)	 20.3/8"
Max Height (cm)	 86	Seat Width (cm)	 39	Max Load (kg)	 150
Minimum Height (cm)	 76	Seat Height from the Ground (cm)	 51	Net Weight (kg)	 6.2