



PIONEER 11

Features

The new streamline sports designed Pioneer range of scooters offer safety, performance & comfort with:

- ✓ Adjustable prestige car like suspension which supplies a smooth ride over varied terrain
- ✓ Built-in speed reduction switch which is activated when turning offering exceptional safety
- ✓ New state of the art anti-glare LCD screen (dash) with introduction of speedo, odometer, clock, temperature gauge, easy to read more accurate bar style battery gauge. The screen also shows when the scooter is in the half speed mode.
- ✓ Dual action parallel tiller adjustment allows the gauges & switches to be visible at all times
- ✓ New look dash layout with easy to see & use switches for indicators, head light, horn, hazards & half speed.
- ✓ Heavy duty 4 pole motor & larger electronic controller that will deliver all the power you need & provide worry free transport
- ✓ Comfortable leather feel contoured high back seat with slide, swivel, recline & armrest height & width adjustment
- ✓ New safety reflectors on each side of the scooter for improved visibility
- ✓ Power saving LED tail lights, brake lights & indicators at the rear
- ✓ Front basket & mirrors standard
- ✓ Unique 12 volt USB fitting (for phone charging)



840DXB-1 PIONEER 11



**New Anti Glare
LCD Screen**



**Adjustable
Suspension**




Handy USB fitting



**Powerful 4 Pole
Motor**

Dimensions



Specifications

Speed up to	10kph
Range up to	35km (dependent on battery size & terrain)
Turn Radius	150cm
Ground Clearance	12cm
Motor	DC 24V, 400W (Rating)/ 1750W (Max)
Controller	Dynamic Rhino 2, 120 amp
Battery	12V/33-50a/h
Charger	5amp off-board
Gradient	8" degrees
Front Wheel	11" Pneumatic
Rear Wheel	11" Pneumatic
Brake	Intelligent, regenerative, electromagnetic
Weight capacity	159kg

*NOTE - 33a/h-50a/h batteries can be fitted

Australian Standards Compliant AS/NZS 3695.2:2013 Pending

Colour Options

- ✓ Red
- ✓ Silver
- ✓ Dark Blue

Peak Care is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.