

Power Pack for manual wheelchair P001F

Features

Easy attachment to manual wheelchairs to allow motorised attendant assistance

Convenient forward & reverse direction switch & speed control, with easy to operate thumb throttle

Portable & lightweight for easy transportation which enables extra mobility anywhere

Twin wheel design which allows much greater traction on grass & gravel



Technical Specs

Motor	280watt
Battery	12V - 20a/h
Charger	3A off board
Drive Wheel	8" PU tyre x 2 pieces
Weight of battery	6.3kg
Weight of base	7.7kg
Total weight	14kg
Weight Capacity	136kg at (6 degree gradient) – 159kg at (0 degree gradient)
Range up to	15km
Fits wheelchair width	16inch – 22inch

Universal mount kit
available as an option