

# PYTHON 4

---

**COLOUR OPTIONS:** MS050PBAU-1 – Blue  
MS050RDAU-1 – Red

- Max range of 48km\*
- Top speed of 6kmph
- Powerful 350W motor
- All round suspension
- Stylish silver wheel hubs
- Delta bars as standard
- Black pneumatic tyres with large 5" rims
- LED front headlight
- Recommended MK batteries: 2 x 12v 40ah\*\*

\* Subject to conditions of use. \*\* Batteries not included.



<b>Length</b>	<b>Height</b>
121cm (47.5")	102cm (40")
<b>Width</b>	<b>Weight</b>
60cm (24")	94kg (207lb))
<b>Max Capacity</b>	<b>Max Gradient</b>
160kg (25st)	8°