

Shoprider AllRounder

ASK FOR THE BRAND SHOPRIDER



NOW WITH
4 POLE
MOTOR

Available Colours



SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER® AllRounder

www.shoprider.com.au *mobility*

The SHOPRIDER™ **AllRounder** is loaded with features that make it a superior choice when selecting a mobility scooter from the Shoprider range. The Shoprider **AllRounder** comes with a powerful 4 pole 1.5hp motor, 330mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold down arm rests along with all the other great features you expect from a Shoprider scooter.

The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake

SPECIFICATIONS

	AllRounder
Length	1300mm
Width	640mm
Height	1200mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
Maximum Load Single User*	159Kgs
Charger	8Amp
Motor Size	1.5Hp (4 pole)
Maximum Speed	10kmph
Battery Capacity (Minimum Amp)	40Amp
Maximum Range #	up to 40km
Wheel Size	330mm
Turning Radius	1440mm
Weight (with battery)	99kg

FEATURES

- Powerful 4 Pole Motor
- Large Captain's Seat
- Adjustable Steering Tiller
- Adjustable Width Arm Rests
- Sliding rotating Seat
- Fold Arm Rests
- Variable Seat height
- Front Shopping Basket
- Horn
- Led Battery Level indicator
- Front Chrome Bumper
- Front & Rear Suspension
- Front & Rear Lights
- Hazard Lights /indicators
- High / Low Speed switch
- Anti Tips Wheels

Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)


WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

AllRounder



SHOPRIDER®
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au