



*For who you are & what you need...*

## Lightweight Rollators



# For who you are & what you need.

At any age, we need to feel good and stay active. Knowing that we can get up and go anywhere, at any time, is vital to your independence and enjoying life to its full potential.

With this in mind, we developed our **Airgo®** rollators. Offering more than just support, they provide convenience, comfort, stability and security. We take great care to create distinct models that offer unique features so that you can choose one that's right for you.

The joy of walking is that it is a safe, easy and effective exercise. An **Airgo®** rollator can help you improve your health and may even increase your energy level over time, as you walk.

There's an **Airgo®** rollator for you. Try one today and feel the difference with every step!



# Our exclusive features make a difference!



## Handles

Only **Airgo®** features **Ergo-form Handgrips** that fit much more comfortably and place your hand at a natural angle for better brake access. They allow you to reach and grip the brakes without shifting, allowing for better control and added safety.

Feather-touch brakes offer easy manoeuvrability.

Available on: eXcursion™, Fusion™, Lightweight Rollator, Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



## Backrests

Our large **Contoured Backrest** provides full back support thanks to its unique profile design. It comfortably encompasses your entire back while providing improved support for your posture.

Available on: eXcursion™, Fusion™ and Vista™.



The **Ergo-Back** is simply the most comfortable rigid rollator backrest available today. It features an ergonomically shaped soft padded cushion made from high density foam.


Available on: Comfort-Plus and Comfort-Plus XWD.



## Seats

The **Ergo-Seat**'s ergonomic shape and high density foam cushion with laminated memory foam on top provides unrivalled seating comfort. Its unique design distributes weight effectively across the entire surface, reducing average pressure measurements by as much as 70% over standard rollator seats.

Available on: Comfort-Plus and Comfort-Plus XWD.



**Ergo-Seat™**  
Pressure map with Piezo resistive technology

**Normal Hard Seat**



## Footrests

Our new **Swing n' Lock Footrests** secure to the frame when folded, allowing for a hassle free experience. Converting from a rollator to transport chair has never been easier!

Available on: Fusion™.



## Basket

Soft Basket (underseat) that does not need to be removed when folding the rollator. With convenient compartments and zipper closure, it is opaque and located under the seat for stability, keeping valuables organized and accessible to the user yet safely out of sight.

Available on: Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



Our covered Soft Basket (side-fold) with reflective stitching doesn't need to be removed when folding your rollator. In addition, it is machine washable and dryable for easy cleaning.

Available on: eXcursion™, Fusion™ and Vista™.



## Wheels

Comfort-grip™ Indoor/Outdoor Wheels are non-marring and offer better grip on carpet along with superior braking action. Solid mags for easy maintenance.

Available on: Comfort-Plus and Comfort-Plus XWD.



Our new Low-Profile Comfort-Grip™ Wheels feature a sport tire design that offers better stability, with low rolling resistance on smooth indoor surfaces. Once outdoors, the outer treads grip in to provide superior control.

Available on: eXcursion™ and Fusion™.



## Brakes

Our new Tool-Free Brake Adjustment makes adjusting your brake cables a breeze. The large knob is senior friendly and FAQ videos are available on our website.

Available on: eXcursion™, Fusion™ and Vista™.

## Height options

We're not all created the same. Our Multi-Height Frames allows comfort and safety to fit even more people by offering our rollators at the seat height that's right for you.

Available on: Adventure 8, Comfort-Plus, Comfort-Plus XWD, Ultra-Light 6, eXcursion™ and Fusion™.



Hemi  
19 ½  
(49 cm)



Standard  
22 ½  
(57 cm)

## Accessories



### Oxygen Bag

703-120

Accommodates size C and M6 oxygen tanks.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus and Comfort-Plus XWD.



### Saddle Bag

703-100

Snaps conveniently onto your rollator's backrest.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



### Side-fold Cane Holder

730-946

Allows you to securely and easily attach a cane to your side-folding rollator.

Fits on: eXcursion™ and Fusion™.



### Cane Holder

730-947

Allows you to securely and easily attach a cane to your rollator.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus and Comfort-Plus XWD.



### Side-fold Cup Holder

730-944

Stay hydrated while walking with this convenient rotatable cup holder.

Fits on: eXcursion™ and Fusion™.



### Apollo LED Safety Light

730-950

Rechargeable LED light featuring 3 usage modes, and 20 Lumens of power.

Fits on: All rollators.



### Side-fold Curb Climber

A01-551

Makes it easier to lift your side-folding rollator over obstacles such as sidewalks.

Fits on: eXcursion and Fusion.





































### Curb Climber

703-050


Makes it easier to lift your rollator over obstacles such as sidewalks.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus, and Comfort-Plus XWD.

# There is an Airgo rollator for your needs!

					
	Colour	Ergo-Form Handgrips	Ergo-Seat and Back	Contoured Backrest	Basket
Lightweight					Wire
Adventure8					Soft closed, underseat
Ultra-Light6					Soft closed, underseat
Comfort-Plus					Soft closed, underseat extra large
Comfort-Plus XWD					
eXcursion X23					Covered, soft, in front of rollator
eXcursion X20	   				
eXcursion X18	 				
eXcursion XWD	 				
eXcursion6					
Fusion F23					Covered, soft, in front of rollator
Fusion F20					
Fusion F18					
Vista	   				Covered, soft, in front



 Seat Adjustment	 Ultra-light Frame	 wheels	 Transport Feature	 MAX. weight capacity	 NET weight
Fixed 22" (56 cm)		Standard 8" (20 cm)		300 lb (136 kg)	16.2 lb (7.3 kg)
Dual height		Standard 8" (20 cm)		300 lb (136 kg)	16.2 lb (7.3 kg)
Multi-height	✓	Standard 6" (15 cm)		250 lb (113 kg)	10.5 lb (7.8 kg)
Dual height		Comfort-Grip 8" (20 cm)		300 lb (136 kg)	17.2 lb (7.8 kg)
				400 lb (181 kg)	17.8 lb (8 kg)
Fixed 23" (58 cm)	✓	Low-Profile Comfort-Grip 8" (20 cm)		300 lb (136 kg)	17.1 lb (7.7 kg)
Fixed 20" (51 cm)				300 lb (136 kg)	15.4 lb (6.9 kg)
Fixed 18" (46 cm)				300 lb (136 kg)	15.4 lb (6.9 kg)
Fixed 20" (51 cm)				400 lb (181 kg)	16.4 lb (7.4 kg)
Fixed 20" (51 cm)		Standard 6" (15 cm)		300 lb (136 kg)	13.5 lb (6 kg)
Fixed 20" (51 cm)	✓	Low-Profile Comfort-Grip 8" (20 cm)	Full 2 in 1	300 lb (136 kg)	17.1 lb (7.7 kg)
Fixed 18" (46 cm)				300 lb (136 kg)	
Fixed 23" (58 cm)				300 lb (136 kg)	18.9 lb (8.6 kg)
Fixed 20" (51 cm)	✓	Low-Profile 8" (20 cm)		300 lb (136 kg)	15 lb (6.8 kg)

# eXcursion<sup>TM</sup> AAC



Ergo-form handgrips



Ultra-light frame



Easy to fold



Recognized by  
the CAOT

## The ultra-light side-fold with superior comfort & stability!

- Full contoured backrest provides superior back support (3).
- Fold by simply pulling the strap (4), the soft covered basket folds with the frame (5).
- A folding lock keeps the unit closed for secure lifting & storing.
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (1).
- Comfort-grip wheels positioned outside the frame for increased stability (6).
- NEW tool-free easy to adjust brake system.
- Weight capacity of 300 lb / XWD weight capacity of 400 lb.

X23 (23" seat)	5'9" - 6'11"	33.5" - 39.5"	23"	17.5"	25"	28.5"
Colour: Pacific Blue (700-917) Net weight (without acc.) 17.1 lb						
X20 (20" seat)	5" - 6'2"	30.5" - 37"	20"	17.5"	25"	28.5"
Colours: Pearl Black (700-920) Pacific Blue (700-921) Cranberry (700-922) Lime (700-923) Net weight (without acc.) 15.4 lb						
X18 (18" seat)	4'3" - 5'5"	29.5" - 35"	18"	17.5"	25"	28.5"
Colours: Pacific Blue (700-926) Pearl Black (700-927) Net weight (without acc.) 15.4 lb						
XWD	5" - 6'2"	30" - 36.5"	20"	19.5"	27"	28.5"
Colours: Pacific Blue (700-928) Pearl Black (700-929) Net weight (without acc.) 16.4 lb						

Patents 7828305, 2587989, 200710103964.9 ADP type 3 MW3GM0065 (eXcursion<sup>TM</sup> X18, X20, X23)  
AADL C225 ADP type 3 AA0000124 (eXcursion<sup>TM</sup> XWD) AADL C064.



# eXcursion<sup>TM</sup> MC



Ergo-form handgrips



Ultra-light frame  
– 6" wheels –



Easy to fold



Recognized by  
the CAOT

Our ultra-light eXcursion frame, now even lighter with 6" wheels!

- Full contoured backrest provides superior back support (3).
- Fold by simply pulling the strap (4), the soft covered basket folds with the frame (5).
- A folding lock keeps the unit closed for secure lifting & storing.
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (1).
- 6" wheels are non-marring with low rolling resistance, and ideal for indoor use. (6).
- NEW tool-free easy to adjust brake system.
- Weight capacity of 300 lb.

X20 (20" seat)	5' - 6'2"	30" - 36"	20"	17.5"	24"	27.5"	
Colours: Pacific Blue (700-919) Net weight (without acc.) 13.5 lb							

Patents 7828305, 2587989, 200710103964.9 • ADP type 2 AA0001355 • AADL C225.

# FUSION<sup>TM</sup> MC



Recognized by  
the CAOT



Comfortable armrests



Easy to fold



Swing n' Lock<sup>TM</sup>  
footrests

Converts effortlessly from a side-folding rollator  
to a transport chair.

- Padded seat and full contoured backrest (1).
- Soft basket doesn't need to be removed when folding (2).
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (3).
- Easy to activate brakes (4).
- Locks in open position and when folded for safer lifting.
- 8" non-marring wheels are ideal for indoor/outdoor use (7).
- Weight capacity: 300 lb (136 kg).
- Seat belt available separately, item 700-863.

Fusion F23	5'9" - 6'11"	33.5" - 39.5"	23"	19.5"	27"	28.5"
Colour: Pearl Black (700-932) Net weight (without acc.) 19.9 lb						
Fusion F20	5' - 6'2"	30" - 36.5"	20"	19.5"	27"	28.5"
Colour: Pacific Blue (700-935) Net weight (without acc.) 17.1 lb						
Fusion F18	4'3" - 5'5"	29.5" - 35"	18"	19.5"	27"	28.5"
Colour: Cranberry (700-938) Net weight (without acc.) 17.1 lb						

Patents 7828305, 2587989, 200710103964.9, 7540527 and patents pending. ADP type 3 AA0000G94.

# Vista™



Our travel-friendly side-folding rollator, with pop-in wheels and a lightweight design.

- Padded seat and adjustable ergo backrest (1).
- Extra-large soft basket with a secure cover, that folds with the frame (2).
- Adjustable ergonomic handles, with integrated parking brakes (3).
- Brakes lock for security when seated (4).
- Easy to fold, and stands when folded. Locks both open and closed for safety.
- Travel friendly pop-in wheels allow it to be stored almost anywhere (5).
- Treaded 8" front wheels for superior steering and rolling comfort (6).

Vista V20	5' - 6'2"	35.5" - 38"	20"	18"	25"	28.5"
Colours: Rocket red (700-906) Royal blue (700-907) Onyx black (700-908) Stellar silver (700-909)						
Net weight (without acc.) 15 lb						

ADP Type 3 AA0001338, AADL C225.

# Ultra-Light6



Ergo-form handgrips



Ultra-light frame



Multi-set height technology

The lightest rollator available:  
a breeze to lift up and take anywhere!

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 6" non-marring wheels are ideal for indoor use (6).
- Reflective trim for added visibility (7).
- Weight capacity: 250 lb (113 kg).

Ultra-Light6	5' - 6'2"	30" - 34"	19 1/4" - 23 1/4"	17"	25"	31"
Colours: Eggplant (700-950) Net weight (without acc.) 10.5 lb						

Patents 7828305, 2587989, 200710103964.9 ADP type 2 MW2GM0050 AADL C225

# Lightweight Rollator



Ergo-form handgrips

An economical choice with all the basic features.

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Removable wire basket (4).
- 8" non-marring wheels are ideal for indoor/outdoor use (5).
- Reflective trim for added visibility (6).
- Weight capacity: 300 lb (136 kg).

Lightweight	5' - 6'2"	33" - 39"	22"	17"	24.5"	27.5"

Colours: Pacific Blue (700-940) Net weight (without acc.) 16.2 lb

Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9 AADL C225.

# Adventure™ 8



Ergo-form handgrips



Ultra-light frame



Dual height frame technology

## Light and dependable!

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 8" non-marring wheels (6).
- Reflective trim for added visibility (7).
- Weight capacity: 300 lb (136 kg).

Standard	5' - 6'2"	33" - 39"	21.5"	17"	24.5"	27.5"
HEMI	4'3" - 6'2"	30" - 36"	18.5"	17"	24.5"	25.5"

Colours: Titanium (700-942) Net weight (without acc.) 16.2 lb

Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9

ADP type 3 MW3GM0060 (Adventure 8) AADL C225.



# COMFORT-PLUS™

## and COMFORT-PLUS™ XWD



Ergo-form handgrips



Dual height  
frame technology



Ergo-Seat & Back



Extra-Large  
shopping basket

Simply the most comfortable rollator on the market.

- Easy-to activate brakes (2).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Comfort-Grip indoor/outdoor 8" wheels (6).
- Reflective trim for added visibility (7).
- Comfort-Plus™ weight capacity: 300 lb (136 kg).
- Comfort-Plus XWDTM weight capacity: 400 lb (181 kg).

Standard	5' - 6'2"	33" - 39"	22.5"	17"	24.5"	27.5"
HEMI	4'3" - 6'2"	30" - 36"	19.5"	17"	24.5"	25.5"
Colour: Iridescent Blue (700-946) Net weight (without acc.) 17.2 lb						
XWD	5' - 6'2"	33" - 39"	22.5"	19"	26.5"	27.5"
	4'3" - 6'2"	30" - 36"	19.5"	19"	26.5"	25.5"
Colour: Iridescent Blue (700-949) Net weight (without acc.) 17.8 lb						

Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9

ADP type 3 MW3GM0050 (Comfort-Plus™) AADL C225.

ADP type 3 MW3GM0055 (Comfort-Plus™ XWD) AADL C064.

# Real Life Engineering

Your Airgo® Rollator is guaranteed to be free from defects in materials and workmanship. This commitment demands that our products are engineered to last, but more importantly, it allows you the confidence and freedom to get up and go, anywhere, anytime!

We also appreciate your feedback so we can learn how to continuously improve our rollators. We call this “Real Life Engineering”.



[www.caot.ca/proudsupporter](http://www.caot.ca/proudsupporter)



[www.airgomobility.com](http://www.airgomobility.com)

Scan to visit our webpage featuring our NEW educational videos.



Proud sponsor of

[www.walker-facts.com](http://www.walker-facts.com)

Please visit to learn more about rollators & walkers.

**Lakeside Mobility**  
**2/1 Metier Linkway**  
**Birtinya QLD 4575**  
**Ph) 5325 1231**

P/N 915-213 rev 11 0516



Drive Medical Canada, Inc.  
256 Aviva Park Drive, Unit #2  
Vaughan, Ontario L4L 9C7

To learn more, call:  
(866) 752-1363 or visit us at

[www.airgomobility.com](http://www.airgomobility.com)