



3000512.D3B02-0125

INVACARE®

# Birdie Compact

*Designed to fold easily and without the need for tools, the Birdie Compact can be stored or transported comfortably when space is limited. The compact base and large castors enables smooth manoeuvring around beds, in bathroom areas and through narrow doorways. With a weight capacity of 150kg, the Birdie Compact is ideal for everyday lifting in the home or aged care facility environment. Compact without compromising on comfort.*

### EASY HANDLING – FOR THE CARER

The Birdie Compact is only 53 cm wide making it the only lifter that is TÜV certified that is this narrow – unique to the market. The ergonomically designed push bar allows the carer to position their hands comfortably when moving the lifter, reducing the likelihood of strain. The 100mm smoothing rolling castors make it easy to push and pull the lifter, reducing the force required to manoeuvre.

### SIMPLE TO FOLD, UNFOLD AND DISMANTLE WITHOUT TOOLS

The Birdie Compact takes up minimal space when stored reducing clutter in confined spaces. When space is limited, the Birdie Compact can even be stored upright.

### TÜV-CERTIFICATION FOR MAXIMUM SECURITY

The Birdie Compact is CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure the Birdie Compact complies with AS 10535.2002 (the Australian standard for hoists) the lifter has been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

### CONFIGURED TO USER NEEDS

Users can enjoy maximum flexibility with a range of spreader bars and a very comfortable sling range, which allows for individual adjustment of the transfer. The new quick release attachment system makes configuration even easier.

### Technical Data

110cm	56.5cm	89cm	54 - 161cm
12.5cm	150kg	35.8kg	80 cycles*
	3 years frame 2 years electronics and motors 1 year battery		

\*Varies depending upon load and stroke.

### Standard:

- 2 Point Spreader Bar
- Manual Leg Adjustment Handle
- Internal Charging

### Accessories:

- 4 Point Spreader Bar
- Pivot Frame
- Weighing Device (250kg SWL)
- Full range of slings



Quick Release Attachments



Tool-less and easy to fold



Store upright to save space