




Standard Deluxe V2

- Recommended Patient Height: Up to 1620mm to 1880mm
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the user
- Adjustable back recline, seat tilt and raising leg rest
- Drop away / Adjustable arm rests and back wings
- Carrflex Extreme multi-stretch covering with sanitised treatment and enhanced durability
- Large castors make the chair easy to manoeuvre
- One directional lock castor to assist with straight line tracking
- Access Vehicle Approved to Australian Standards
- Custom colours available

Air comfort

		MAX USER WEIGHT (Kg)	OVERALL LENGTH (MM)	OVERALL HEIGHT (MM)	PRODUCT WEIGHT (Kg)	SEAT WIDTH (MM)	SEAT DEPTH (MM)	FLOOR TO SEAT HEIGHT (MM)	OVERALL WIDTH (MM)
Standard Deluxe V2 -Day Chair	AC59111	180	1130	1260	45	470	550	570	680
Standard Deluxe V2 - Day Chair - Electric	AC59112	180	1130	1260	45	470	550	570	680



Small Deluxe V2

- Recommended Patient Height: Up to 1420mm.
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the user.
- Adjustable back recline, seat tilt and raising leg rest.
- Drop away / Adjustable arm rests and back wings.
- Carrflex Extreme multi-stretch covering with sanitised treatment and enhanced durability.
- Large castors make the chair easy to manoeuvre.
- One directional lock castor to assist with straight line tracking.
- Access Vehicle Approved to Australian Standards
- Custom colours available

Air comfort

		MAX USER WEIGHT (Kg)	OVERALL LENGTH (MM)	OVERALL HEIGHT (MM)	PRODUCT WEIGHT (Kg)	SEAT WIDTH (MM)	SEAT DEPTH (MM)	FLOOR TO SEAT HEIGHT (MM)	OVERALL WIDTH (MM)
Small Deluxe V2 - Day Chair	AC59114	180	1000	1130	40	430	400	550	670