

Karma Ergo Lite Deluxe

Code: KM-2512F14



TAC
Approved



DVA
Approved



Australian
Standards



Features

- Frame:** Ultra light weight aluminium folding
- Footrest:** Swing away removable
Height Adjustable
- Armrest:** Ergonomic Pads
- Wheels:** Fixed
14" x 1.75" spoke
Solid Tyres
- Castors:** Solid 6"
- Brakes:** Attendant system
- Seat:** Fire retardant upholstery
- Back:** Fire retardant upholstery
folding
- Warranty:** 2 Years - Frame



S-shaped ergonomic seating system

The Ergo Lite Deluxe is the lightest transit wheelchair in our product range. The wheelchair has Karma's patented S-Ergo seating system which promotes pressure relief and helps prevent the user from sliding forward.

The Ergo Lite Deluxe has a double cross brace for strength, durability and quality you would expect from Karma Wheelchairs.

Options



Anti Tip Bars



Travel Bag



Rear Carry Bag



Seat Belt



Handle Extension

Specifications

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
Ergo Lite 16"x16"	410mm	420mm	480mm	420mm	100kg	8.7kg
Ergo Lite 18"x16"	450mm	420mm	480mm	420mm	100kg	8.8kg

Karma
Wheelchairs