

# GelBodies

Gelbodies are sold in pairs

## Heel & Ankle Protection

	Product Code	Dimensions
Small	GB030	08-10" / 20-25cm
Medium	GB031	10-12" / 25-30cm
Large	GB032	12-14" / 30-36cm

Protects both the heel and ankle from the development of pressure sores.

Measure around the ankle above the maleoli (the two protuberances). Choose the larger size if measurement is between sizes.



## Heel & Elbow Protection

	Product Code	Dimensions
Small	GB010	08-10" / 20-25cm
Medium	GB011	10-12" / 25-30cm
Large	GB012	12-14" / 30-36cm
X-Large	GB013	14-16" / 36-41cm

Protects both the heel and elbow from the development of pressure sores.

Measure around the calf or upper arm. Choose the larger size if measurement is between sizes.



## Shin Protection

	Product Code	Dimensions
Medium	GB021	10-13" / 25-33cm
Large	GB022	13-16" / 33-41cm

Protects against skin breakage in the shin area.

Measure around the calf. Choose the larger size if measurement is between sizes.



## Knee Protection

	Product Code	Dimensions
Small	GB040	09-11" / 23-28cm
Medium	GB041	11-13" / 28-33cm
Large	GB042	13-16" / 33-41cm

Protect the knee from the development of pressure sores.

Measure above the knee. Choose the larger size if measurement is between sizes.



Cleaning Instructions:

Coolband: See washing symbols.

Gel Pad: Remove from coolband. Wipe with a damp cloth using standard cleaning agent. It is not advised to soak gel over prolonged periods.

