

# Karma Eagle

Code: KM-8020QD24



TAC  
Approved



DVA  
Approved



Australian  
Standards



## Features

- Frame:** Light weight aluminium folding
- Footrest:** Height Adjustable  
Swing away removable
- Armrest:** Swing Back  
Ergonomic Pads  
Height Adjustable
- Wheels:** Quick release  
24" x 1 3/8" spoke  
Pneumatic Tyres
- Castors:** Solid 8"
- Brakes:** Pull to lock
- Seat:** Height adjustable  
Fire retardant upholstery
- Back:** Removable padding  
Fire retardant upholstery  
Height adjustable
- Warranty:** 10 Years - Frame

Transit Approved

The Karma Eagle is a lightweight, durable and reliable manual wheelchair that makes it the right choice for both users and carers.

The Eagle offers an impressive range of standard features with an extensive selection of personal customising available. The Eagle wheelchair is easy to fold and transport and is the logical choice for people seeking a quality mobility product at an affordable price.

The use of aircraft-grade materials and unique design of the cross-bars and link bars provide this model with outstanding strength and durability.

## Specifications

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
Eagle 14"x16"	360mm	410mm	500mm	400mm	130kg	14.2kg
Eagle 16"x16"	410mm	410mm	500mm	400mm	130kg	14.4kg
Eagle 16"x18"	410mm	460mm	500mm	400mm	130kg	14.5kg
Eagle 17"x16"	430mm	410mm	500mm	400mm	130kg	14.5kg
Eagle 18"x16"	460mm	410mm	500mm	400mm	130kg	14.6kg
Eagle 18"x18"	460mm	460mm	500mm	400mm	130kg	14.6kg
Eagle HD 19"x18"	490mm	460mm	520mm	400mm	160kg	15.5kg
Eagle HD 20"x18"	510mm	460mm	520mm	400mm	160kg	15.7kg
Eagle HD 22"x18"	550mm	460mm	520mm	400mm	160kg	15.9kg