

Karma S-Ergo 125

Code: KM-S125Q



**TAC
Approved**



**DVA
Approved**



**Australian
Standards**

Features



- Frame:** Light weight aluminium
Folding
- Footrest:** Height Adjustable
Swing away removable
- Armrest:** Swing back
Ergonomic Pads
- Wheels:** Quick release
24" x 1 3/8" spoke
Solid Tyres
- Castors:** Solid 7"
- Brakes:** Push to lock
Attendant system
- Seat:** Removable padding
Fire retardant upholstery
- Back:** Removable padding
Fire retardant upholstery
Folding
- Warranty:** 5 Year - Frame



**S-shaped ergonomic
seating system**

Karma's S-Ergo 125 is equipped with the Unique S-Ergo shaped seat. The S-Ergo seat provides pressure relief and prevents the user from sliding forward. In addition to Karma's patented S-Ergo seating system the Ergo 125 also has seat upholstery and pads treated with Aegis Microbe shield, to reduce stains and odour. The stream-line oval tubes not only make the wheelchair look great but also add to the strength of the wheelchair structure, which is incredibly light. Made from aircraft aluminium the frame breaks down to just over 10kg.

Options



Travel Bag



Anti Tip Bars



Elevating Legrests



Oxygen Holder



Folding Headrest



Stump Support



Handle Extension

Specifications

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
S-Ergo 125 16" x 17"	410mm	430mm	500mm	400mm	115kg	14.2kg
S-Ergo 125 18" x 17"	450mm	430mm	500mm	400mm	115kg	14.5kg
S-Ergo 125 20" x 17"	510mm	430mm	500mm	400mm	115kg	14.8kg