

Karma Transit II

Code: KM-8520F12



TAC
Approved



DVA
Approved



Australian
Standards



Features

- Frame:** Light weight aluminium folding
- Footrest:** Height Adjustable
Swing away removable
- Armrest:** Swing Back
Ergonomic Pads
Height Adjustable
- Wheels:** Fixed
12" x 2 1/4" mag
Solid Tyres
- Castors:** Solid 8"
- Brakes:** Push to lock
Attendant system
- Seat:** Height adjustable
Fire retardant upholstery
- Back:** Fire retardant upholstery
Height adjustable
folding
- Warranty:** 10 Years - Frame

Transit Approved

The Karma Transit II is a lightweight, durable and reliable manual wheelchair that makes it the right choice for both users and carers.

The Transit II is the attendant propelled version of the Karma Eagle and offers an impressive range of standard feature with an extensive selection of personal customising available.

The use of aircraft-grade materials and unique design of the cross bars and link bars provide this model with outstanding strength and durability.

Specifications

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
Transit II 16"x16"	410mm	410mm	520mm	400mm	130kg	13.5kg
Transit II 18"x16"	460mm	410mm	520mm	400mm	130kg	14.3kg
Transit II 18"x18"	460mm	460mm	520mm	400mm	130kg	14.3kg
Transit II 20"x18"	510mm	460mm	520mm	400mm	160kg	15.2kg
Transit II 22"x18"	550mm	460mm	520mm	400mm	160kg	15.7kg