

Seating / Day Chairs

Chair Electric Lift and Recline Chair



Mobile Compact Lift Chair

- Total surface air pressure system
- Removable seat and back cushions
- Electric lift and recline mechanism
- 100mm front & rear castors with rear locking castor
- Wall Saver mechanism
- Wooden hand knuckles
- Chargeable battery
- Proven pressure relief for long-term seating
- Assisted sitting, lifting, and reclining
- Easy to operate and manoeuvre the client while they are seated in the chair
- Easy to clean, repair and maintain
- Custom colours available

Air comfort

		MAX USER WEIGHT (KG)	OVERALL HEIGHT (MM)	SEAT WIDTH (MM)	SEAT DEPTH (MM)	FLOOR TO SEAT HEIGHT (MM)	BACKREST HEIGHT (MM)	OVERALL WIDTH (MM)
Mobile Compact Lift Chair Small - Single Motor	AC59098	140	1060	480	440	480	680	750
Mobile Compact Lift Chair Medium - Single Motor	AC59087	140	1060	500	510	520	740	750
Mobile Compact Lift Chair Large - Single Motor	AC59099	140	1060	600	510	520	800	820
Mobile Compact Lift Chair Small - Twin Motor	AC59100	140	1060	480	440	520	680	750
Mobile Compact Lift Chair Medium - Twin Motor	AC59090	140	1060	500	510	520	740	750
Mobile Compact Lift Chair Large - Twin Motor	AC59101	140	1060	600	510	520	800	820