

Ocean VIP

COMFORTABLE, ERGONOMIC AND STABLE



The Aquatec Ocean series of height adjustable shower commode chairs offers the right model for all demands.



Yes, you can.®

AQUATEC



Ocean VIP

- Stable, stainless steel powder coated frame
- Seat tilt can be individually adjusted (tilt angle 0° - 35°)
- Stainless steel gas spring
- Extremely stable even at maximum tilt angle
- Seat height adjustable in 3 steps, no tools required
- Very comfortable, ergonomically contoured seat
- Machine-washable adjustable backrest with quick-release buckles – can be easily adjusted to the user

The Ocean Vip is available as follows:

- Standard version with castors
- Standard version with castors and toilet pan holder
- Standard version with castors, toilet pan holder and soft seat
- XL version (more space between the armrests) with castors
- The Ocean Vip can be retrofitted with XL armrests

Features



Seat Tilt angle 0° - 35°.



Height and depth-adjustable headrests included.



Easy to operate brakes on 3 castors, 1 castor with directional lock to help push the chairs down straight floors.



Seat height adjustable in 3 steps, no tools required.

Tech Specs

Ocean VIP	Seat height 51 - 61cm	Seat width 50cm	Total width 59.5cm	Total depth 106cm	Total height without headrest 108.5 - 118.5cm	Max. user weight 150kg	Unit weight 23.5kg



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com

Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



Yes, you can.®