

Peak Coolsion Foam & Gel Cushion

Product Code: 9346376030809

Features:

- High Quality PU Foam
- Aquous Non Toxic Gel
- For Wheelchair & Standard Seating
- Gel Pack Sandwiched and Sealed
- Helps Skin Remain Cool Dry
- Removable Cover
- Easy Maintenance



| Product Specifications | |
|------------------------|-------|
| Length | 45cm |
| Width | 40cm |
| Height | 7.5cm |
| Weight | 2.3kg |
| ILD (IFD 40%) | 120N |
| Tensile Strength | 70KPA |

| Product Specifications | |
|------------------------|----------------------------|
| Foam Type | PU Foam |
| Cover Type | PU Fabric & Anti Slip Base |
| Density | 25kg/m ³ |
| Tear Strength | 70KPA |
| Resistance | 40% |
| Flammability | CA117 |