

Karma S-Ergo 125 Transit

Code: KM-S125T



**TAC
Approved**



**DVA
Approved**



**Australian
Standards**



Features

Frame:	Light weight aluminium Folding
Footrest:	Height Adjustable Swing away removable
Armrest:	Swingback Ergonomic Pads
Wheels:	Fixed 14" x 1.75" mag Solid Tyres
Castors:	Solid 7"
Brakes:	Push to lock Attendant system
Seat:	Removable padding Fire retardant upholstery
Back:	Removable padding Fire retardant upholstery Folding
Warranty:	5 Year -Frame

**S-shaped ergonomic
seating system**

Karma's S-Ergo 125 is equipped with the Unique S-Ergo shaped seat. The S-Ergo seat provides pressure relief and prevents the user from sliding forward. In addition to Karma's patented S-Ergo seating system the Ergo 125 also has seat upholstery and pads treated with Aegis Microbe shield, to reduce stains and odour. The stream-line oval tubes not only make the wheelchair look great but also add to the strength of the wheelchair structure, which is incredibly light.

Options



Travel Bag



Anti Tip Bars



Elevating Legrests



Oxygen Holder



Folding Headrest



Stump Support



Handle Extension

Specifications

Description	Seat Width	Seat Depth	Seat Height	Back Height	Max User Weight	Chair Weight
S-Ergo 125 Transit 16" x 17"	410mm	430mm	510mm	400mm	115kg	13.1kg
S-Ergo 125 Transit 18" x 17"	450mm	430mm	510mm	400mm	115kg	13.1kg
S-Ergo 125 Transit 20" x 17"	510mm	430mm	510mm	400mm	115kg	13.3kg