

Shoprider VENTURER

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®
www.shoprider.com.au

mobility

VENTURER

The SHOPRIDER™ VENTURER sets new standards for design and function for the SHOPRIDER™ range. This scooter offers a modern streamline design with many updated features including an easy to use tiller dash, full LED lighting with super bright headlight, gas adjustable tiller and 265 mm pneumatic tyres. The **Venturer** is also fitted with comfort ride suspension and a high back seat. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by a powerful 1.0hp motor, rear wheel drive and electromagnetic brake.

SPECIFICATIONS

	VENTURER
Length	1250mm
Width	600mm
Height	1160mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	8Amp
Motor Size	1.0Hp
Maximum Speed	9.5Kph
Battery Capacity (Minimum Amp)	50Amp
Maximum Range #	up to 40km
Wheel Size	265mm
Turning Radius	1440mm
Weight (with battery)	95kg

FEATURES

High Back Seat
 Gas Adjustable Steering Tiller
 Easy to Use Dash Panel
 LED Lighting
 Super Bright Headlight
 Variable Seat Height
 265mm Pneumatic Tyres
 Sliding & Rotating Seat
 Adjustable Width Arm Rests
 Front & Rear Suspension
 Front Shopping Basket
 Hazard Lights / indicators
 Horn
 Anti Tips Wheels
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)


WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

VENTURER



SHOPRIDER®
www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
 Or visit our Web Site at www.shoprider.com.au