

Peak Neptune Self-Propelled Wheelchair 16x16

Product Code: 9346376060721 (Red)













Features:

- Red Aluminum Frame
- Drop-Back Handle
- Flip-Up Armrest
- Detachable Footrest
- Durable PU Rear Wheel
- Easily Folds



ARTG Listing no.

- 313946

| | | | | | |
|-------------------|--|----------------------------------|--|---------------------------|--|
| Length (cm) |  100 | Backrest Height (cm) |  42 | Front Wheel Diameter (cm) |  20/8" |
| Width (cm) |  56.5 | Seat Depth (cm) |  40 | Rear Wheel Diameter (cm) |  60/23" |
| Height (cm) |  92 | Seat Width (cm) |  41 | Max Load (kg) |  100 |
| Folded Width (cm) |  28 | Seat Height from the Ground (cm) |  47 | Net Weight (kg) |  13.1 |