

## Peak Taima Rollator with Gutter Arms - Small

Product Code: C4507-B-FX

## **Features:**

- Au Std Testing No. 492166
- One Hand Folding Mechanism
- Secure Brake Design
- Mesh Bag Included
- Cane Holder Included
- Foldable

## **ARTG Listing:**

227601



Length (cm)		63	Folded Width (cm)		22	Front Castor Diameter (cm) 20.3/8"
Width (cm)		58	Seat Depth (cm)		22	Rear Castor Diameter (cm) 20.3/8"
Max Heigh (cm)	† <u>†</u>	130	Seat Width (cm)	<del> </del>	39	Max Load (kg) KG 150
Minimum Height (cm)		114	Seat Height from the Ground(cm)		51	Net Weight (kg) KG 9