

Modular Lift Recliner



ModularLift Recliner

The Australian Made Modular lift recliner can be adjusted on site, giving the user a customised chair in a matter of minutes.



- 3 motors Independent backrest, independent footrest and tilt-in-space recline
- · Lift & tilt or vertical lift configuration options
- · Optional electric headrest
- Home button to return the chair to starting position
- · 2 headrest style options

- · Storage pocket on both sides of the chair
- Australian made mechanism
- Steel frame
- · Adjustable On Site
- · Available in wide range of Versotela fabrics
- 100% Manufactured in Australia by Oscar Furniture

SPECIFICATIONS













	Seat Height	Seat Width	Seat Depth	Arm Height	Backrest Height	Footrest Length
STARTING MEASUREMENT	445mm	500mm	450mm	95mm	760mm	445mm
ADJUSTMENT	+90mm	+80mm	+100mm	+135mm	+70mm	+40mm
MAXIMUM MEASUREMENT	535mm	580mm	550mm	230mm	830mm	485mm

MODIFICATIONS

SEAT DEPTH

ARM HEIGHT

FOAM

BACKREST HEIGHT

HEAD AND ARM REST COVERS

SEAT WIDTH

SEAT HEIGHT

ADJUSTABLE FOOTREST

VERTICAL LIFT

FABRIC SELECTION



All VersoTela fabrics come standard with a suite of high-end specifications, such as, antimicrobial coating, mold resistance, scratch resistance and fire resistance.

Our wide range of collections each boast unique properties that are carefully designed to suit every customer's needs. Whether it be built for use in a care facility, or as a functional but stylish addition to the living room, our VersoTela fabric options are sure to tick all the boxes when it comes to crafting your perfect chair.

Please get in touch with one of our friendly sales representatives for fabric samples and to find out more about which fabric is right for you.









Contact us today: sales@oscarfurniture.com.au 03 5381 1404

6 King Drive, Horsham VIC 3400

www.oscarfurniture.com.au

